Bacterial Viral Filter for Pulmonary Function Testing Equipment

English

Compact Bacterial/Viral Filter with electrostatic filter media, nose clip and mouthpiece. Filters are available in various equipment end-fitting adaptable to various pulmonary function testing (PFT) equipment only.

Trade Name: NEUMOFILT® Ergo

Models:

NEUMOFILT® Ergo	Basic UDI-DI: 843655723NEUMOFILT3G
NEUMOFILT® Ergo Pack Comfort	Basic UDI-DI: 843655723NEUMOFILTPCRR
NEUMOFILT® Ergo Pack Plus	Basic UDI-DI: 843655723NEUMOFILTPPSK
NEUMOFILT® Ergo Pack Comfort Plus	Basic UDI-DI: 843655723NEUMOFILTPCPVQ

References for model:

MODEL NEUMOFILT® Ergo

I.D: Internal diameter O.D: Outer diameter	REFERENCE	UDI- DI
I.D. 30mm O.D. 32 mm	1532	8436557230003
I.D. 45.5 mm O.D. 48 mm	1533	8436557230164
I.D. 34 mm O.D. 37.5 mm	1534	8436557230324
I.D. 27.6mm O.D. 30.4 mm	1535	8436557230485
I.D. 28.6mm O.D. 30.9 mm	1536	8436557230645

MODEL NEUMOFILT® Ergo Pack Plus

I.D: Internal diameter O.D: Outer diameter	REFERENCE	UDI- DI
I.D. 30mm O.D. 32 mm	1532 N03	8436557231031
I.D. 45.5 mm O.D. 48 mm	1533 N03	8436557230201
I.D. 34 mm O.D. 37.5 mm	1534 N03	8436557231086
I.D. 27.6mm O.D. 30.4 mm	1535 N03	8436557231161
I.D. 28.6mm O.D. 30.9 mm	1536 N03	8436557231284

MODEL NEUMOFILT® Ergo Pack Comfort

I.D: Internal diameter O.D: Outer diameter	REFERENCE	UDI- DI
I.D. 30mm O.D. 32 mm	1532 M	8436557230089
I.D. 45.5 mm O.D. 48 mm	1533 M	8436557230249
I.D. 34 mm O.D. 37.5 mm	1534 M	8436557230409
I.D. 27.6mm O.D. 30.4 mm	1535 M	8436557230560
I.D. 28.6mm O.D. 30.9 mm	1536 M	8436557230720

MODEL NEUMOFILT® Ergo Pack Comfort Plus

I.D: Internal diameter O.D: Outer diameter	REFERENCE	UDI- DI
I.D. 30mm O.D. 32 mm	1532 N03M	8436557231079
I.D. 45.5 mm O.D. 48 mm	1533 N03M	8436557231246
I.D. 34 mm O.D. 37.5 mm	1534 N03M	8436557230447
I.D. 27.6mm O.D. 30.4 mm	1535 N03M	8436557230607
I.D. 28.6mm O.D. 30.9 mm	1536 N03M	8436557231208

INTENDED USE:

Device is intended for the prevention of cross-contamination during pulmonary function testing (PFT): spirometry, plethysmography and Helium Dilution Measurement).

INTENDED USERS

Healthcare professionals with specific training and instruction on respiratory function testing. NO specific training for the device use is required.

INTENDED PATIENT POPULATION

Patients undergoing respiratory function tests. There are no restrictions beyond the discretion of the health professional to submit the patient to these tests. The health professional must consider the contraindications associated with respiratory function tests.

ENVIRONMENT OF USE: Hospital, Sub-acute facilities, Clinics, Physician offices.

COMPATIBILITY: the device is a general device intended for use in combination with equipment for respiratory function testing such as spirometers or plethysmographs. There are no general restrictions for use with pulmonary function testing (PFT) equipment, but there are restrictions regarding diameter compatibility between such PFT equipment, and the diameters offered by the device. Connection is made by applying pressure.

INDICATIONS FOR USE:

Performance of respiratory function tests. For use with pulmonary function testing (PFT) equipment only, to filter air between the patient's exhaled air and the testing equipment. Filter and components are disposable, single patient, single PFT procedure.

The device is designed to be used in pulmonary function test equipment such as spirometry, plethysmography and helium dilution measurement. For the proper use of device follow the PFT manufacturer's equipment's indications.

The acceptability of a particular PFT procedure is the responsibility of the clinical staff. Special attention should be given when testing elderly patients, children, pregnant women and people with disabilities and regarding to warnings, indications and contraindications of PFT procedure.

User Instructions Bacterial/Viral Filter for pulmonary function testing equipment:

1. Connect the cylindrical part of the filter to the pulmonary function testing equipment, (see Fig.1)
2. The patient must perform the test by breathing (inhaling and exhaling) through the ergonomic side of the filter, (see Fig.1.)
3. Maximal utilization time required time for carrying out the following tests: spirometry, plethysmography, Lung Volumes by dilution.
4. SINGLE USE ONLY FOR SINGLE PATIENT PROCEDURE. NEVER REUSE THE FILTER EVEN WITH THE SAME PATIENT. RE-USE THE FILTER CAN CAUSE CROSS INFECTION BETWEEN PATIENTS AS WELL AS INADEQUATE FILTRATION.

User Instructions Noseclip (when applicable)

1. Remove from package.
2. Test prior to use by opening and closing the nose clip at least once.
3. Attach to the patient's nose.
4. SINGLE USE ONLY FOR SINGLE PATIENT PROCEDURE. NEVER REUSE THE NOSE CLIP EVEN WITH THE SAME PATIENT. RE-USE THE NOSE CLIP CAN CAUSE AN INADEQUATE OCCLUSION OF THE NOSTRILS.

User Instructions Rubber mouthpiece: (when applicable).

1. Remove from package.
2. Place and press strongly the mouthpiece over the filter patient side (see Fig 1 and 2)
3. SINGLE USE ONLY FOR SINGLE PATIENT PROCEDURE. NEVER REUSE THE MOUTH PIECE EVEN WITH THE SAME PATIENT. RE-USE. THE RUBBER MOUTHPIECE CAN CAUSE PRODUCT DEFORMATION.

Warnings:

1. In case of any serious incident that has occurred in relation to the device, it should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established. For more report information see: https://health.ec.europa.eu/medical-devices-sector/new-regulations/contacts_en.
2. Verify that the diameter of the filter connection matches the diameter of the pulmonary function testing equipment.
3. Store filters in a clean, dry environment away from sources of heat or direct light. Handle the device with clean hands to avoid contamination. Users are recommended to wear disposable gloves when handling the filter.
4. Before use, verify that the device is not damaged, perforated or contaminated. Discard filters with any signs of damage or deterioration.
5. Incorrectly placing the filter in the pulmonary function testing equipment may result in inaccurate readings or lack of adequate protection against contamination
6. The duration of regular PFT are around 15-20 minutes.
7. Dispose of the device at the end of its useful life in accordance with local regulations on medical device waste. Do not waste it in general waste, as it may represent a risk to health and the environment.

Precautions:

1. The filter, nose clip and mouthpiece are disposable intended for one patient during single pulmonary function testing. They should not be reused even for the same patient in different pulmonary function tests, as this increases the risk of cross-contamination and may compromise the filter's effectiveness.
2. The filter neither the nose clip and mouthpiece cannot be cleaned or sterilized after use. Device should be discarded after the PFT procedure.

Contraindications:

1. Potential relative and absolute contraindications for the pulmonary function testing (PFT) should be evaluated carefully by healthcare professional for particular patient.
2. Patients with pacemakers do not have an absolute contraindication to lung function testing. However, it is important to inform the physician about the presence of the pacemaker before performing the lung function test.
3. Not for use in patients known to be allergic to any of the device components.

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English:

The IFU shall not be provided in electronic format. Single copy of the IFU shall be provided in paper format into the box in compliance with GSPR 23.1(e). Where multiple devices are supplied to a single user and/or location, a single copy of the instructions for use may be provided if so, agreed by the purchaser who in any case may request further copies to be provided free of charge.

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